


Guaiacol-D3

Chemical Properties

CAS No. :	74495-69-5
Formula:	C7H8O2
Molecular Weight:	127.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Guaiacol-D3 is a deuterated form of Guaiacol (T1301), which is a phenolic compound known for its ability to inhibit COX-2 expression and NF-κB activation in response to LPS stimulation, exhibiting anti-inflammatory properties.
Targets(IC50)	Others,NF-κB,COX

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.8641 mL	39.3205 mL	78.6411 mL
5 mM	1.5728 mL	7.8641 mL	15.7282 mL
10 mM	0.7864 mL	3.9321 mL	7.8641 mL
50 mM	0.1573 mL	0.7864 mL	1.5728 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481