

YZ-268

Chemical Properties

CAS No. :


Formula: C₄₉H₅₇N₉O₇

Molecular Weight: 884.03

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	YZ-268 is a PROTAC that targets HDAC.
Targets(IC50)	PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1312 mL	5.6559 mL	11.3118 mL
5 mM	0.2262 mL	1.1312 mL	2.2624 mL
10 mM	0.1131 mL	0.5656 mL	1.1312 mL
50 mM	0.0226 mL	0.1131 mL	0.2262 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481