

Dth

## Chemical Properties


CAS No. :

Formula: C30H29F2N5O6

Molecular Weight: 593.58

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dth is a proteolysis-targeting chimera (PROTAC) composed of the fluorescent RNA probe DFHBI and the immunomodulatory compound thalidomide.
Targets(IC50)	NF-κB,PROTACs,Early 2 Factor (E2F)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6847 mL	8.4235 mL	16.8469 mL
5 mM	0.3369 mL	1.6847 mL	3.3694 mL
10 mM	0.1685 mL	0.8423 mL	1.6847 mL
50 mM	0.0337 mL	0.1685 mL	0.3369 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481