

P22D

Chemical Properties

CAS No. :


Formula: C₄₇H₅₄ClN₁₃O₉S

Molecular Weight: 1012.53

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | P22D is a PROTAC that targets BCR-ABL. |
| Targets(IC50) | PROTACs |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.9876 mL | 4.9381 mL | 9.8763 mL |
| 5 mM | 0.1975 mL | 0.9876 mL | 1.9753 mL |
| 10 mM | 0.0988 mL | 0.4938 mL | 0.9876 mL |
| 50 mM | 0.0198 mL | 0.0988 mL | 0.1975 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481