

SIAIS126


Chemical Properties

CAS No. :

Formula: C44H47ClFN9O7

Molecular Weight: 868.35

Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	SIAIS126 is a PROTAC targeting EGFR.
Targets(IC50)	PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1516 mL	5.758 mL	11.5161 mL
5 mM	0.2303 mL	1.1516 mL	2.3032 mL
10 mM	0.1152 mL	0.5758 mL	1.1516 mL
50 mM	0.023 mL	0.1152 mL	0.2303 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481