

SPB5208

Chemical Properties

CAS No. :


Formula: C44H47N7O9

Molecular Weight: 817.89

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|------------------------------------------------------------------------|
| Description | SPB5208 is a PROTAC designed to target BTK (Bruton's tyrosine kinase). |
| Targets(IC50) | PROTACs |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.2227 mL | 6.1133 mL | 12.2266 mL |
| 5 mM | 0.2445 mL | 1.2227 mL | 2.4453 mL |
| 10 mM | 0.1223 mL | 0.6113 mL | 1.2227 mL |
| 50 mM | 0.0245 mL | 0.1223 mL | 0.2445 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481