

UI-EP002

## Chemical Properties

CAS No. :


Formula: C60H87N5O16S

Molecular Weight: 1166.42

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	UI-EP002 is a PROTAC compound designed to target ER (estrogen receptor).
Targets(IC50)	PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8573 mL	4.2866 mL	8.5732 mL
5 mM	0.1715 mL	0.8573 mL	1.7146 mL
10 mM	0.0857 mL	0.4287 mL	0.8573 mL
50 mM	0.0171 mL	0.0857 mL	0.1715 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481