

BSJ-05-002

Chemical Properties

CAS No. :


Formula: C51H59F3N8O8

Molecular Weight: 969.06

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|------------------------------------|
| Description | BSJ-05-002 serves as a FAK PROTAC. |
| Targets(IC50) | PROTACs |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.0319 mL | 5.1596 mL | 10.3193 mL |
| 5 mM | 0.2064 mL | 1.0319 mL | 2.0639 mL |
| 10 mM | 0.1032 mL | 0.516 mL | 1.0319 mL |
| 50 mM | 0.0206 mL | 0.1032 mL | 0.2064 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481