

TJA-107

Chemical Properties

CAS No. :


Formula: C₄₃H₄₆ClN₇O₆

Molecular Weight: 792.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TJA-107 is an AR PROTAC.
Targets(IC50)	PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2621 mL	6.3106 mL	12.6212 mL
5 mM	0.2524 mL	1.2621 mL	2.5242 mL
10 mM	0.1262 mL	0.6311 mL	1.2621 mL
50 mM	0.0252 mL	0.1262 mL	0.2524 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481