

## 4-Methyl-1-phenyl-2-pentanol

## Chemical Properties

CAS No. : 7779-78-4


Formula: C<sub>12</sub>H<sub>18</sub>O

Molecular Weight: 178.27

Keep away from moisture

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Methyl-1-phenyl-2-pentanol is an alcohol containing a benzene ring that is used in organic synthesis and biochemical experiments.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6095 mL	28.0473 mL	56.0947 mL
5 mM	1.1219 mL	5.6095 mL	11.2189 mL
10 mM	0.5609 mL	2.8047 mL	5.6095 mL
50 mM	0.1122 mL	0.5609 mL	1.1219 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481