

## Thieno[2,3-b]thiophene

## Chemical Properties

CAS No. : 250-84-0


Formula: C<sub>6</sub>H<sub>4</sub>S<sub>2</sub>

Molecular Weight: 140.23

Storage: Keep away from direct sunlight, Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thieno[2,3-b]thiophene is a thiophene that can be used in organic optoelectronic materials and biochemical experiments.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.1311 mL	35.6557 mL	71.3114 mL
5 mM	1.4262 mL	7.1311 mL	14.2623 mL
10 mM	0.7131 mL	3.5656 mL	7.1311 mL
50 mM	0.1426 mL	0.7131 mL	1.4262 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481