

## p-Tolyl Phenylacetate

## Chemical Properties

CAS No. : 101-94-0


Formula: C<sub>15</sub>H<sub>14</sub>O<sub>2</sub>

Molecular Weight: 226.27

Keep away from moisture

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	p-Tolyl Phenylacetate is an organic ester compound commonly used in flavor preparation, organic synthesis and biochemical experiments.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4195 mL	22.0975 mL	44.195 mL
5 mM	0.8839 mL	4.4195 mL	8.839 mL
10 mM	0.4419 mL	2.2097 mL	4.4195 mL
50 mM	0.0884 mL	0.4419 mL	0.8839 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Api AM, et al. RIFM fragrance ingredient safety assessment, p-tolyl phenylacetate, CAS Registry Number 101-94-0. Food Chem Toxicol. 2025 Apr 3:115423.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481