

## CCL4L1 Protein, Human, Recombinant

### General Information

Synonyms:	LAG1;SCYA4L1;CCL4L1;Macrophage inflammatory protein 1-beta (MIP-1-beta);SCYA4L2; Monocyte adherence-induced protein 5-alpha;Lymphocyte activation gene 1 protein (LAG-1); Small-inducible cytokine A4-like;C-C motif chemokine 4-like;CCL4L2;CCL4L
Protein Construction:	24-92 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q8NHW4
Molecular Weight:	7.8 kDa (Predicted)
AA Sequence:	APMGSDPPTACCFSTARKLPRNFVVDYYETSSLCSQPAVVFQTKRGKQVCADPSESWVQEYVVDLELN

### QC Testing

Biological Activity:	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using human CCR5 transfected murine BaF3 cells is less than 2.0 ng/ml, corresponding to a specific activity of $> 5.0 \times 10^5$ IU/mg.
Purity:	$>97\%$ as determined by SDS-PAGE.
Endotoxin:	$< 1.0$ EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 $\mu$ m filtered PBS, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at  $-80^{\circ}\text{C}$ . For reconstituted protein solutions, the solution can be stored at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481