

Lymphotactin/XCL1 Protein, Human, Recombinant

General Information

| | |
|-----------------------|---|
| Synonyms: | Lymphotaxin;SCYC1;LTN;Cytokine SCM-1;XC chemokine ligand 1;ATAC;Lymphotactin;Small-inducible cytokine C1;C motif chemokine 1;XCL1;SCM-1-alpha |
| Protein Construction: | 23-114 aa |
| Species: | Human |
| Expression Host: | E. coli |
| Accession: | P47992 |
| Molecular Weight: | 10.2 kDa (Predicted) |
| AA Sequence: | GSEVSDKRTCVSLTTQRLPVSRIKTYTITEGSLRAVIFITKRLKVCADPQATWVRDVVRSMDRKSNTNRNMIQ TKPTGTQSTNTAVTLTG |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human T-lymphocytes is in a concentration of 10-100 ng/ml. |
| Purity: | >98% as determined by SDS-PAGE. |
| Endotoxin: | < 1.0 EU/µg of the protein as determined by the LAL method. |
| Formulation: | Lyophilized from a 0.2 µm filtered PBS, pH 7.4 |

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481