

## Mitofusin-2 Protein, Human, Recombinant (His)

## General Information

Synonyms: KIAA0214;CPRP1;MFN2;Mitofusin-2;Transmembrane GTPase MFN2

Protein Construction: 1-757 aa

Species: Human

Expression Host: in vitro E. coli expression system

Accession: O95140

Molecular Weight: 87.9 kDa (Predicted)

MSLLFSRCNSIVTVKKNKRHMAEVDNASPLKHFVTAKKKINGIFEQLGAYIQESATFLEDTYRNAELDPVTTEEQ  
VLDVKGYSKVRGISEVLARRHMKVAFFGRTSNGKSTVINAMLWDKVLPSPGIGHTTNCFLRVEGTDGHEAFLL  
TEGSEEKRSKTVNQLAHALHQDKQLHAGSLVSMWPNSKCPDLLKDDLVLMDSPGIDVTTELDSWIDKFCLD  
ADVFVLVANSESTLMQTEKHFFHKVSRNIFILNNRWDASASEPEYMEVRRQHMERCSTSLVDELGVV  
DRSQAGDRIFVSAKEVLNARIQKAQGMPEGGALAEQFQVRMFEFQNFERRFEECISQSAVKTKFEQHTVR  
AA Sequence: AKQIAEAVRLIMDSLHMAAREQQVYCEEMREERQDRLKFIDKQLELLAQDYKLRKQITEEVERQVSTAMAEI  
RRLSVLVDDYQMDFHPSVVLKVYKNELHRHIEEGLGRNMSDRCSITNSLQTMQQDMIDGLKPLLPVSVR  
SQIDMLVPRQCFSLNYDLNCDKLCADFQEDIEFHFSLGWTMLVNRFLGPKNSRRALMGYNDQVQRPIPLTPA  
NPSMPPLPQGSLTQEEFMVSMVTGLASLTSRTSMGILVGGVWVWKAQVWRLIALSFGLYGLLYVYERLTWTTK  
AKERAFKRQFVEHASEKLQLVISYTGNSCSHQVQQELSGTFAHLCQQVDVTRENLEQEIAMNKKIEVLDSLQ  
SKAKLLRNKAGWLDSELNMFTHQYLQPSR

## QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Purity: >85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.

Formulation: If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 34 Washington Street, Wellesley Hills, MA 02481