

## NT-4 Protein, Human, Recombinant (Active)

## General Information

Synonyms:	NTF5;Neurotrophin-5 (NT-5);Neurotrophic factor 4;NTF4;Neurotrophin-4;NT-4
Protein Construction:	81-210 aa
Species:	Human
Expression Host:	E. coli
Accession:	P34130
Molecular Weight:	14.1 kDa (Predicted)
AA Sequence:	M+GVSETAPASRRGELAVCDVSGWVTDRRTAVDLRGREVEVLGEVPAAGGSPLRQYFFETRCKADNAEEG GPGAGGGGCRGVDRRHVVSECKAKQSYVRALTADAQGRVGVWRWIRIDTACVCTLLSRTGRA

## QC Testing

Biological Activity:	Fully biologically active when compared to standard. The ED50 as determined by the dose-dependent induction of choline acetyl transferase activity in rat basal forebrain primary septal cell cultures is less than 50 ng/ml, corresponding to a specific activity of $> 2.0 \times 10^4$ IU/mg
Purity:	$>97\%$ as determined by SDS-PAGE.
Endotoxin:	$< 1.0$ EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 $\mu$ m filtered PBS, pH 5.5

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at  $-80^{\circ}\text{C}$ . For reconstituted protein solutions, the solution can be stored at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481