

## NOV/CCN3 Protein, Human, Recombinant

## General Information

Synonyms:	CCN family member 3;IGFBP9;NOVH;NOV;Cellular communication network factor 3;CCN3;Nephro blastoma-overexpressed gene protein homolog;Insulin-like growth factor-binding protein 9 (IBP-9;IGF-binding protein 9;IGFBP-9);Protein NOV homolog (NovH)
Protein Construction:	28-357 aa
Species:	Human
Expression Host:	E. coli
Accession:	P48745
Molecular Weight:	36.2 kDa (Predicted)
AA Sequence:	M+QVAATQRCPQCPGRCPATPPTCAPGVRAVLGDCSCCLVCARQRGESCSDLEPCDESSGLYCDRSADPS NQTGICTAVEGDNCVFDGVIYRSGEKFPQSCFKQCTCRDGGIGCVPRQLDVLPEPNCAPRKEVVPGECC KWICGPDEEDSLGGLTLAAYRPEATLGVEVSDSSVNCIEQTTEWTACSKSCGMGFSTRVTNRNRQCEMLKQT RLCMVRPCEQEPEQPTDKKGGKCLRTRKSLKAIHLQFKNCTSLHTYKPRFCGVCSDGRCCTPHNTKTIQAEFQ CSPGQIVKKPVMVIGTCTCHTNCNPKNNEAFLQELELKTTRGKM

## QC Testing

Biological Activity:	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using murine Balb/c 3T3 cells is less than 1.0 µg/ml, corresponding to a specific activity of > 1000 IU/mg.
Purity:	>95% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.6, 150 mM NaCl.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481