

PTPRT Protein, Human, Recombinant (hFc)

General Information

| | |
|-----------------------|---|
| Synonyms: | Receptor-type tyrosine-protein phosphatase rho (RPTP-rho); Receptor-type tyrosine-protein phosphatase T; KIAA0283; PTPRT; R-PTP-T |
| Protein Construction: | 26-747 aa |
| Species: | Human |
| Expression Host: | HEK293 Cells |
| Accession: | O14522 |
| Molecular Weight: | 109.2 kDa (Predicted) |
| AA Sequence: | QSAAGGCSFDEHYSNCGYSVALGTNGFTWEQINTWEKPMLDQAVPTGSFMMVNSSGRASGQKAHLLLP KENDTHCIDFHYYFSSDRSSPGALNVYVKVNGGPQGNPVWNVSGVVTEGWVKAELAISTFWPHFYQVIFES VSLKGHPGYIAVDEVRVLAHPCRKAPHFLRLQNEVNVGQNFQCIAGGKWSQHDKLWLQWNGRDTAL MVTRVNVHRRFSATVSVADTAQRSVSKYRCVIRSDGGSGVSNYAELIVKEPPTPIAPPELLAVGATYLVWIKPNA NSIIGDGPILKEVEYRTTGTWAETHIVDSPNYKLWHLDPDVEYEIRVLLTRPGEGGTGPPGPPLTTRTKCADP VHGPQNVEIVDIRARQLTLQWEPFGYAVTRCHSYNLTVQYQYVFNQQQYEAEEVIQTSSHYTLRGLRPFMTIR LRLLSNPEGRMESEELVVQTEEDVPGAVPLESIQGGPFEEKIYIQWKPNETNGVITLYEINYKAVGSLDPSAD LSSQRGKVFKLNRNETHHLFVGLYPGTTYSFTIKASTAKGFGPPVTTRIATKISAPSMPEYDTPNETDITITVM LKPAQSRGAPVSVYQLVVKEERLQKSRAADIIECFVSPVSYRNASSLDSLHYFAAELKPANLPVTQPFTVGD NKTYNGYWNPLSPLKSYSIYFQALS KANGETKINCVRLATKGASTQNSNTVEPEKQVDNTVKM |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Activity has not been tested. It is theoretically active, but we cannot guarantee it. |
| Purity: | >85% as determined by SDS-PAGE. |
| Endotoxin: | < 1.0 EU/μg of the protein as determined by the LAL method. |
| Formulation: | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481