

SVEP1 Protein, Human, Recombinant (His)

General Information

Synonyms:	CCP module-containing protein 22; Selectin-like osteoblast-derived protein (SEL-OB); von Willebrand factor type A; CCP22; C9orf13; SELOB; EGF and pentraxin domain-containing protein 1; Sushi, von Willebrand factor type A, EGF and pentraxin domain-containing protein 1; Serologically defined breast cancer antigen NY-BR-38; SVEP1; Polydom
Protein Construction:	736-827 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q4LDE5
Molecular Weight:	17.6 kDa (Predicted)
AA Sequence:	GDFICTPDNTGVNCTLTCLLEGYDFTEGSTDKYYCAYEDGVWKPTYTTEWPDCAKKRFANHGFKSFEMFYKAA RCDDTDLMKKFSEAFETTLG

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	>95% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481