

## TAF12 Protein, Human, Recombinant (His & SUMO)

### General Information

Synonyms:	TAFII20;TAF15;TAF2J;Transcription initiation factor TFIID subunit 12;TAF12;Transcription initiation factor TFIID 20/15 kDa subunits (TAFII-20/TAFII-15;TAFII20/TAFII15)
Protein Construction:	1-161 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q16514
Molecular Weight:	33.9 kDa (Predicted)
AA Sequence:	MNQFGPSALINLSNFSSIKPEPASTPPQGSMANSTAVVKIPGTPGAGGRLSPENNQVLTKKKLQDLVREVDP NEQLDEDVEEMLLQIADDFIESVVTAAACQLARHRKSSTLEVKDVQLHLERQWNMWIPGFGSEEIRPYKKA CTT EAHKQRMALIRKTTKK

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	>90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

### Preparation and Storage

Reconstitution:	Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Stability & Storage:	Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481