

## Annexin V/ANXA5 Protein, Mouse, Recombinant (His &amp; Myc)

## General Information

Synonyms:	Annexin-5;Vascular anticoagulant-alpha (VAC-alpha);Thromboplastin inhibitor;Anx5; Endonexin II;Annexin A5;Annexin V;Placental anticoagulant protein 4 (PP4);Anxa5; Calphobindin I (CPB-I);Anchorin CII;Lipocortin V;Placental anticoagulant protein I (PAP-I)
Protein Construction:	2-319 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	P48036
Molecular Weight:	43.2 kDa (Predicted)
AA Sequence:	ATRGTVTDFPGFDGRADAEVLRKAMKGLGTDEDSILNLLTSRSNAQRQEIAQEFKTLFGRDLVDDLKSELTKG FEKLIVAMMKPSRLYDAYELKHALKGAGTDEKVLTEIIASRTPEELSAIKQVYEEYGSNLEDDVVGDTSGYYQ RMLVLLQANRDPDTAIDDAQVELDAQALFQAGELKWTDEEKFITIFGTRSVSHLRRVFDKYM TISGFQIET IDRETSGNLEQLLLAVVKSIRSIPAYLAETLYYAMKGAGTDDHTLIRVVVSRSEIDLFNIRKEFRKNFATSLYSMIK GDTSGDYKKALLLLCGGEDD

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	>95% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481