

## CCL9 Protein, Mouse, Recombinant (Active)

### General Information

Synonyms:	Mrp2;C-C motif chemokine 9;Macrophage inflammatory protein 1-gamma (MIP-1-gamma); Scya9;Scya10;Ccl9;Small-inducible cytokine A9;Macrophage inflammatory protein-related protein 2 (MRP-2);CCF18
Protein Construction:	22-122 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	P51670
Molecular Weight:	11.6 kDa (Predicted)
AA Sequence:	QITHATETKEVQSSLKAQQGLEIEMFHMGGFQDSSDCCLSYNSRIQCSRFIGYFPTSGGCTRPGLIIFISKRGFQVC ANPSDRRVQRCIERLEQNSQPRTYKQ

### QC Testing

Biological Activity:	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human neutrophils is in a concentration range of 0.1-10 ng/ml.
Purity:	>95% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 µm filtered PBS, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481