

## NTNG1 Protein, Mouse, Recombinant (His & Myc)

### General Information

Synonyms: Kiaa0976; Netrin-G1; Lmnt1; Laminet-1; Ntn1

Protein Construction: 19-510 aa

Species: Mouse

Expression Host: E. coli

Accession: Q8R4G0

Molecular Weight: 63.2 kDa (Predicted)

AA Sequence: VMQPYLFWGHYDVCKSLIYTEEGKVWDYTACQPESTDMTKYLKVKLDPPDITCGDPPESFCAMGNPYMCN  
NECDASTPELAHPPELMFDFEGRHPSTFWQSATWKEYPKPLQVNITLSWSKTIELTDNIVITFESGRPDQMILE  
KSLDYGRTWQPYQYYATDCLHAFHMDPKSVKDLQHTVLEICTEYSTGYSTNSKIIHFEIKDRFAFFAGPRLR  
NMSLYGQLDTTKLRDFFVTDLRIRLLRPAVGEIFVDELHLARYFYAISDIKVRGRCKCNLHATSCLYDNSKL  
TCECEHNTTGPDCGKCKKNYQGRPWSYLPKGTANTCIPSISSIGNCECFGHSNRCSYIDLLNTVICVSK  
HNTRGQHCELRLGYFRNASALDDENVCIECYCNPLGSIHDCRNGSGFCECKTGTGPKCDECLPGNSWYY  
GCQPNVCDNELLHCQNGGTCQNNVRCACPDAYTGILCEKLRCEEAGSCGS

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Purity: >85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.

Formulation: If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.  
If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

### Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481