

ALDH1A2 Protein, Mouse, Recombinant (His)

General Information

Synonyms:	RaLDH2;RALDH 2;Aldh1a7;Retinaldehyde-specific dehydrogenase type 2 (RALDH(II)); Aldehyde dehydrogenase family 1 member A2 (ALDH1A2);Retinal dehydrogenase 2;Aldh1a2
Protein Construction:	1-518 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q62148
Molecular Weight:	62.6 kDa (Predicted)
AA Sequence:	MTSSEIAMPGEVKADPAALMASLQLLPSTPNLEIKYTKIFINNEWQNSESGRVFPVCNPATGEQVCEVQEAD KVDIDKAVQAARLAFSLGSVWRRMDASERGRLLDKLADLVERDRATLATMESLNGGKPFLLQAFYIDLQGVIKT LRYAGWADKIHGMTIPVDGDYFTTRHEPIGVCGQIIPWNFLLMFTWKIAPALCCGNTVVIKPAEQTPLSAL YMGALIKEAGFPVGVNLPYGGPTAGAAIASHIGIDKIAFTGSTEVGKLIQEAGRSNLKRVTLLEGGKSPNIIF ADADLDYAVEQAHQGVFFNQGCCTAGSRIFVEESIYEEFVKRSVERAKRRIVGSPFDPTEQGPQIDKKQYN KVLELIQSGVAEGAKLECGGKGLGRKGGFFIEPTVFSNVTDDMRIAKEEIFGVPVQEILRFKTMDEVIERANNSDFG LVAAVFTNDINKALMVSSAMQAGTVWINCYNALNAQSPFGGFKMSGNGREMGEFGLREYSEVKTIVTKIPQ KNS

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	>85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:	Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Stability & Storage:	Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>

Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.
-----------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481