

Thrombopoietin Protein, Mouse, Recombinant

General Information

Synonyms:	MGDF;C-mpl ligand;Thpo
Protein Construction:	22-195 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	P40226
Molecular Weight:	18.7 kDa (Predicted)
AA Sequence:	SPVAPACDPR LLNLLRDSH LLHSRLSQCP DVDPLSIPVL LPAVDFSLGE WKTQTEQSKA QDILGAVSLL LEGVMAARGQ LEPSCLSLL GQLSGQVRL LGALQGLLGT QLPLQGRRTA HKDPNALFLS LQQLLRGKVR FLLLVEGPTL CVRRTLPTTA VPSSTSQLLT LNKF

QC Testing

Biological Activity:	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using human MO7e cells is less than 1.0 ng/ml, corresponding to a specific activity of $> 1.0 \times 10^6$ IU/mg.
Purity:	$> 95\%$ as determined by SDS-PAGE and HPLC.
Endotoxin:	< 0.1 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 μ m filtered concentrated solution in PBS, pH 7.4, 5% trehalose.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water containing 0.1% BSA. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C . For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481