

## Anti-Mouse PD-1 Antibody (S-5001)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-MousePD-1Antibody (S-5001) is an IgG1 antibody inhibitor that targets mouse PD-1.
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481