

Anti-Mouse LPAM-1/Integrin $\alpha 4\beta 7$ Antibody (DATK32)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse LPAM-1/Integrin $\alpha 4\beta 7$ Antibody (DATK32) is a rat-derived IgG2a, κ antibody inhibitor targeting mouse LPAM-1/Integrin $\alpha 4\beta 7$.
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481