

Anti-ADAM9 Antibody (MGA021)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-ADAM9 Antibody (MGA021) (IMGC-936 antibody) is a humanized antibody expressed in CHO cells, targeting ADAM9. It features hulgG1 heavy chains and huk light chains, with an estimated molecular weight (MW) of 147.62 kDa. For its isotype control, refer to HumanIgG1kappa, Isotype Control.
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481