

NGM-120

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 47.77 kDa

Storage: Store at low temperature,  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NGM-120 is a human-derived antibody that targets GFRAL and can be used in studies of receptor-antibody interactions involving the GFRAL signaling pathway.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0209 mL	0.1047 mL	0.2093 mL
5 mM	0.0042 mL	0.0209 mL	0.0419 mL
10 mM	0.0021 mL	0.0105 mL	0.0209 mL
50 mM	0.0004 mL	0.0021 mL	0.0042 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481