

NGM-120

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NGM-120 is a human-derived antibody expressed in CHO cells and targets GFRAL. It features a huIgG1 heavy chain and a huk light chain, with an estimated molecular weight (MW) of 145.7 kDa. The isotype control for NGM-120 is HumanIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481