

PY159

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	PY159 is a humanized antibody expressed in CHO cells, targeting TREM1/CD354. It features a human IgG1 heavy chain and a human κ light chain, with a predicted molecular weight (MW) of 143.66 kDa. For its isotype control, see Human IgG1 kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481