

## Sapelizumab

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Sapelizumab is a humanised antibody targeting IL-6Ra/CD126, containing huIgG2-type heavy chains and huk-type light chains, used as an IL-6 antagonist for neuromyelitis optica spectrum disorders (NMOSD).
Targets(IC50)	IL Receptor,Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481