

Rezorstobart

Chemical Properties

CAS No. : 2762201-86-3

Formula:

Molecular Weight:

Storage: Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rezorstobart is a humanised monoclonal antibody targeting KLRB1 (Killer cell lectin-like receptor subfamily B, member 1), applicable for immunological research.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481