

Solrikitung

Chemical Properties

CAS No. : 2768721-24-8

Formula:

Molecular Weight:

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Solrikitung is a humanized IgG1κ monoclonal antibody targeting CRLF2, with its humanized IgG1 kappa Isotype Control serving as the isotype control counterpart.
Targets(IC50)	Interleukin,JAK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481