

Anti-Mouse CXCR3/CD183 Antibody (CXCR3-173)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-Mouse CXCR3/CD183 Antibody (CXCR3-173) is an IgG antibody inhibitor derived from a hamster originating from the United States, specifically targeting mouse CXCR3/CD183.
Targets(IC50)	CXCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481