

Anti-Mouse CXCL9/MIG Antibody (MIG-2F5.5)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse CXCL9/MIG Antibody (MIG-2F5.5) is a monoclonal antibody targeting CXCL9/MIG. Anti-Mouse CXCL9/MIG Antibody (MIG-2F5.5) can be used in immunological research and studies on cancers such as prostate cancer.
Targets(IC50)	CXCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481