

Anti-Mouse RANKL/CD254 Antibody (IK22/5)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Anti-MouseRANKL/CD254Antibody (IK22/5) is a rat-derived IgG2a, κ antibody inhibitor targeting the mouse RANKL/CD254. |
| Targets(IC50) | RANKL/RANK |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481