

Anti-Mouse/Human/Rat/Monkey/Hamster/Canine TGF β Antibody (1D11.16.8)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse/Human/Rat/Monkey/Hamster/Canine/Bovine TGF- β Antibody (1D11.16.8) is an IgG antibody inhibitor targeting TGF- β in multiple species including mouse, human, rabbit, monkey, hamster, dog, and cow. For information on the isotype control, refer to Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	TGF-beta/Smad

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481