

## Anti-Mouse IL-7R $\alpha$ /CD127 Antibody (A7R34)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse IL-7R $\alpha$ /CD127 Antibody (A7R34) is an IgG2a, $\kappa$ antibody inhibitor derived from rats that targets mouse IL-7R $\alpha$ /CD127.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481