

Anti-Mouse/Rat/Human PD-L1 Antibody (368A.4H1)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse/Rat/HumanPD-L1Antibody (368A.4H1) is a murine-derived IgG1, κ antibody inhibitor targeting PD-L1 in mice, rats, and humans.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481