

## Anti-Mouse FAP Antibody (73.3)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	The compound Anti-MouseFAPAntibody (73.3) is an IgG1, κ antibody inhibitor derived from mice, designed to target and inhibit mouse FAP.
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481