

## Anti-Mouse/Rat/Bovine VLDL-R Antibody (IgG-6A6)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse/Rat/Bovine VLDL-R Antibody (IgG-6A6) is a murine-derived IgG1, $\lambda$ antibody inhibitor targeting VLDL-R in mice, rats, and bovines.
Targets(IC50)	LDLR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481