

## Anti-Mouse MHC Class I (H-2Kd, H-2Dd) Antibody (34-1-2S)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse MHC Class I (H-2Kd, H-2Dd) Antibody (34-1-2S) is an IgG2a, κ antibody inhibitor derived from mice and targets murine MHC Class I.
Targets(IC50)	MHC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481