

Anti-Mouse/Human KLRG-1 Antibody (2F1)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse/Human KLRG-1 Antibody (2F1) is an IgG antibody inhibitor originating from Syrian hamster, targeting mouse or human KLRG-1.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481