

## Anti-Mouse IL-12 p40 Antibody (C17.8)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse IL-12p40 Antibody (C17.8) is a rat-derived IgG2a, κ inhibitor targeting mouse IL-12p40.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481