

## Anti-Mouse/Human IL-7 Antibody (M25)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse/Human IL-7 Antibody (M25) is a mouse-derived IgG2b antibody inhibitor that targets mouse or human IL-7.
Targets(IC50)	IL Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481