

Anti-Mouse IL-15 Antibody (AIO.3)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse IL-15 Antibody (AIO.3) is a rat-derived IgG2a, λ antibody inhibitor targeting mouse IL-15.
Targets(IC50)	IL Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481