

NGM-313

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	NGM-313 is a CHO-expressed human antibody that targets Klotho Beta. It features a hulgG1 heavy chain and a huk light chain, with a predicted molecular weight (MW) of 150 kDa. For an isotype control of NGM-313, refer to HumanIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481