

OS-2966

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	OS-2966 is a humanised de-immunised monoclonal antibody that targets and inhibits the adhesion receptor subunit Integrin b1/ITGB1/CD29, and can be used for the treatment of malignant gliomas, glioblastomas, and astrocytomas.
Targets(IC50)	Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481