

## Anti-CYR61/CCN1 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature,  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-CYR61/CCN1 Antibody is a humanized antibody that targets CYR61/CCN1. It is suitable for studying mechanisms such as extracellular matrix signaling, angiogenesis-related pathways, cell adhesion, and the biology of the tumor microenvironment.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481